




2 Witney Way
 Boldon Business Park
 Tyne and Wear
 NE35 9PE
 England

New Sterile Disposable Alger Brush Burrs

The following products are now available from
 Altomed Limited.




To place an order contact Customer Services on
 Tel: (0) 191 519 0111 Fax: (0) 191 519 0283



NEW STERILE DISPOSABLE ALGERBRUSH BURRS

During metal working procedures such as welding, sanding or grinding, it is possible that a welding spark / piece of metal can be ejected from the machine into the patient's cornea. If this occurs it is essential that all traces of the foreign body and resulting rust are removed. The Algerbrush has been specially designed to remove rust spots from a patient's cornea.

We are pleased to announce the availability of CE marked single use Altomed burrs for the popular Algerbrush 'rust ring' removal device



Manufactured by Alger – Reusable

A4955 Algerbrush
 Rust Ring Remover
 Machine

A4958 Replacement
 Chuck



Background

Due to concerns about cleaning reusable burrs and also the emergence of BSE and vCJD as an issue in the hospital environment, ophthalmic professionals have expressed the wish to be supplied with disposable burrs.

Product

The new single use burr is a premium quality 0.5mm diameter burr which is identical in specification to the reusable Alger burr. It is available in boxes of 5 or 25.

The burrs are rigorously quality checked prior to packing and sterilization by gamma irradiation.

Manufactured by Altomed – Official Single Use Burrs

A4959 0.5mm Sterile Algerbrush Burrs - Pack of 5
A4960 0.5mm Sterile Algerbrush Burrs - Pack of 25

A4959 1.0mm Sterile Algerbrush Burrs - Pack of 5
A4960 1.0mm Sterile Algerbrush Burrs - Pack of 25



*Altomed also offer a
Repair Service
 for broken Algerbrush Rust Ring
 Removers*