



2 Witney Way  
 Boldon Business Park  
 Tyne and Wear  
 NE35 9PE  
 England




*Altomed presents*

# QuickSet

## LSO Silicone Oil with Quickset, the convenient and clean delivery kit

The following products are now available from Altomed Limited.

To place an order contact Customer Services on  
 Tel: (0) 191 519 0111 Fax: (0) 191 519 0283



Handling and dispensing silicone oil during vitreo retinal surgery has, in the past, created problems relating to;

- Messy silicone oil residue on gloves and instruments -
- Wastage of oil when decanting into a syringe -
- Air bubbles -
- Time taken in decanting setup & dispensing -
- Cost of separate components -

The LSO Silicone Oil and QuickSet 'Combo Pack' eliminates many of these issues by giving you one sterile pack that contains a pre-filled syringe, connected to a pressure tube kit that will fit onto most automated oil pump machines

### 'Combo Pack'

Use the Combo Pack and

- Save time in set-up
- No pouring of oil required
- No messy oil residue from decanting
- No wasted oil left in a vial
- Fits most machines

### Ordering information

LC1000 Combo Pack with 10ml 1000c oil  
 LC5000 Combo Pack with 10ml 5000c oil



The Oil

## LSO™ SILICONE OIL

Ultrapure LSO Silicone Oil from Labtician, Canada.  
 World leader in retinal detachment implants



LSO Silicone Oil in a syringe

'Combo Pack'  
 LSO Silicone Oil in a syringe with  
 QuickSet interface tubing