



2 Witney Way
 Boldon Business Park
 Tyne and Wear
 NE35 9PE
 England

The Altomed Microwash Q-Rinse Machine And Accessories.

The following products are available from
 Altomed Limited.

To place an order contact Customer Services on
 Tel: (0) 191 519 0111 Fax: (0) 191 519 0283

The *Q-Rinse*
 System Is Fully



Marked

A11070 – Microwash Q-Rinse System

Includes:

- Q-Rinse Machine
- A11072 - Bottle Assembly and Tubing.
- A11074 - Power Adapter and Lead
- A11076 - Wire Brush
- A11046 - Assorted Adapters pack of 8
- Instruction Manual



Q-Rinse Accessories.

All fully interchangeable with the AOI Q-Rinse Machine



A11072
Bottle Assembly
Complete



A11044A
Bottle Only



A11044B
Metal Cap



A11044C
4 x Cap Gaskets




A11044D
Full Cap
Assembly



A11044E
Full Tubing Set



A11044F
Bottle Connector
"O" Rings



A11042
4 x Interface
Tubes and Luers



A11074
Power Adapter
and Lead



A11076
Wire and Nylon
Brush



A11046
8 x Assorted
Adapters



A11047C
Elbow Connector



A11049
10 x Male to
Female Luers



A11049A
10 x Female to
Male Luers

Full List of Q-Rinse Devices and Spare Parts.

Replacements:

<u>Code</u>	<u>Description</u>
A11072	Q-Rinse Bottle Assembly – Complete
A11074	Q-Rinse Power Adapter and Lead
A11076	Q-Rinse Wire Brush

Bottle Parts:

A11042	Q-Rinse Spare Interface Tubing
A11044A	Q-Rinse Bottle Only – No Cap or Tubing
A11044B	Q-Rinse Bottle Cap Only – No Tubing
A11044C	Q-Rinse Bottle Cap Gaskets – Pack of 4
A11044D	Q-Rinse Bottle Cap Complete – Includes Tubing
A11044E	Q-Rinse Tubing Set – Front and Back
A11044F	Q-Rinse Bottle Connector “O” Rings

Adaptors:

A11046	Q-Rinse Assorted Adapters – Pack of 8
A11047	Q-Rinse Dual Rinse Tubing “Y” Fitting
A11047C	Q-Rinse Air Input Tube Elbow Connector
A11046	Q-Rinse Assorted Adapters – Pack of 8
A11047	Q-Rinse Dual Rinse Tubing “Y” Fitting
A11047C	Q-Rinse Air Input Tube Elbow Connector

Machine Parts:

A11040C	Q-Rinse Circuit Board
A11043P	Q-Rinse Pump and Motor Unit
A11045	Q-Rinse Case – Upper
A11045L	Q-Rinse Case – Lower
A11047G	Q-Rinse Water (Green) Cycle Switch
A11047S	Q-Rinse Air (Blue) Stop Switch
A11047VA	Q-Rinse Pinch Valve – Air
A11047VF	Q-Rinse Pinch Valve – Fluid
A11049	Q-Rinse Luer Adapters – Female to Female
A11049A	Q-Rinse Luer Adapters – Male to Male

In the event you require Machine Parts, please contact Altomed Limited if you require any details on their installation.