

DECLARATION OF MRI COMPATIBILITY

Labtician Ophthalmics, Inc's Silicone Retinal Implants (solid silicone and sponge) are manufactured from 100% long term implantable grade silicone (polydimethylsiloxane) and as such are MRI compatible.

Labtician Tantalum Clips used in conjunction with Labtician Silicone Retinal Implants are manufactured from pure tantalum metal (a non-ferrous metal) and as such these are also MRI compatible.

The applicable units are:

S 1981-4	Oval Sponge	S 1986-3.5	Oblong Sponge	S 3004	Silicone Tire 9mm
S 1981-5	Oval Sponge	S 1986-4	Oblong Sponge	S 3010	Silicone Tire 6mm
S 1981-7	Oval Sponge	S 1987-7	Accessory Sponge	S 3012	Boat
S 1981-7.5	Oval Sponge	S 1987-8	Accessory Sponge	S 3013	Silicone Strip
S 1982-1.5	Round Sponge 1.5mm	S 2950	Circling Band	S 3014	Silicone Tire 7mm
S 1982-2	Round Sponge 2mm	S 2965	Silicone Strip	S 3014 L	Silicone Tire 7mm WG
S 1982-2.5	Round Sponge 2.5mm	S 2967	Silicone Strip	S 3016	Silicone Tire 10mm
S 1982-3	Round Sponge 3mm	S 2968	Silicone Strip	S 3017	Silicone Tire
S 1982-4	Round Sponge 4mm	S 2969	Circling Band	S 3018	Silicone Sleeve
S 1982-5	Round Sponge 5mm	S 2970	Silicone Strip	S 3019	Silicone Sleeve
S 1982-5T	Round Tunnel Sponge 5mm	S 2971	Silicone Strip	S 3020	Meridional Implant
S 1982-7	Round Sponge 7mm	S 2976	Silicone Button	S 3021	Wedge
S 1982-7.5	Round Sponge 7.5mm x 50mm	S 2977	Silicone Tire 10mm	S 3022	Wedge
S 1983-2.3	Grooved Sponge	S 2978	Meridional Implant - Flat	S 3025	Wong Meridional Implant
S 1983-2.5	Grooved Part Thickness Sponge	S 2986	Silicone Tire 7mm	S 3071	Silicone Sleeve
S 1983-2.8	Grooved Part Thickness Sponge	S 2987	Circling Band	S 3083	Labtician Oval Sleeve™
S 1983-3	Grooved Sponge	S 2988	Meridional Implant	S 3084	Labtician Oval Sleeve™
S 1983-3.2	Grooved Sponge	S 2989	Boat	S 3093	Labtician Buckling Sleeve™
S 1983-3.5	Grooved Sponge	S 2991	Silicone Tire 9mm	S 3094	Labtician Buckling Sleeve™
S 1983-4	Grooved Sponge	S 2992	Silicone Strip	S 4050	Silicone Strip
S 1983-5	Grooved Sponge	S 2993	Pad	S 4055	Silicone Strip
S 1983-7.5	Grooved Sponge	S 2994	Silicone Tire 6.0mm	S 4060	Silicone Strip
S 1984-2.5	Part Thickness Sponge	S 2995	Silicone Tire 8.5mm	S 4065	Silicone Strip
S 1984-3	Partial Thickness Sponge	S 2996	Silicone Tire 10mm	S 5020	Silicone Lace™ 2mm
S 1984-3.3	Partial Thickness Sponge	S 2996 L	Silicone Tire 10mm LG	S 5025	Silicone Lace™ 2.5mm
S 1984-5	Partial Thickness Sponge	S 2997	Silicone Tire 12.5mm	S 566	Double Tantalum Clip 2.0mm
S 1984-7.5	Partial Thickness Sponge	S 2998	Silicone Strip	S 567	Double Tantalum Clip 2.5mm
S 1985-5	Scholda Partial Thickness	S 2999	Silicone Tire 7.0mm		
S 1985-7	Scholda Partial Thickness	S 2999 L	Silicone Tire Wide Groove		

Phil Cuscuna
Director of QA/RA

Date