



2 Witney Way  
Baldon Business Park  
Tyne and Wear  
NE35 9PE  
England

# Reusable Bipolar Cables and Connectors

The following products are available from  
Altomed Limited.

All Devices  
Are



To place an order contact Customer Services on  
Tel: (0) 191 519 0111 Fax: (0) 191 519 0283

**British Standard (BS) #1**  
Fits old Downs, Siemens, Erbe



Pictures  
not to  
scale

**British Standard (BS) #2 (EU Block)**  
Accepts block end Eschman, Aesculap, etc.



**Half Round**  
Fits old Eschmann, G.U.



**Two Pin (American Fitting)**  
Accepts twin round pins, Codman etc.



**Two Plug**  
Twin fixed plugs at either 22mm or 28.6mm apart



**Two Flat Pin**  
Accepts block end Eschmann, Aesculap, etc.



## Assorted Reusable Altomed Limited Bipolar Cables

- A5075: Millennium (BS#1) to Stortz
- A5076: Millennium (BS#1) to 2 Flat Pin
- A5078: Half Round to Stortz Fitting
- A5080-22: 2 Plugs fixed at 22mm apart to 2 Pin
- A5080-28.6: 2 Plugs fixed at 28.6mm apart to 2 Pin
- A5081: 2 Pin to Half Round
- A5082: BS#1 to 2 Pin
- A5084: BS#1 to BS#2
- A5085: Half Round to BS#2
- A5086-22: 2 Plugs fixed at 22mm apart to BS#2
- A5086-28.6: 2 Plugs fixed at 28.6mm apart to BS#2
- A5090-22: 2 Plugs fixed at 22mm apart 2 Flat Pin
- A5090-28.6: 2 Plugs fixed at 28.6mm apart to 2 Flat Pin
- A5095: Stortz Fitting to 2 Pin
- A5096: Stortz Fitting to BS#2
- A5098: BS to 2 Pin
- 2058-71: Oertli Blue Ring to BS
- 2058-72: Oertli Blue Ring to 2 Pin
- 2063-71: Oertli Red Ring to BS
- 2063-72: Oertli Red Ring to 2 Pin

Contact Altomed Limited Customer Services for other cable specifications available including **Single Use Options**.

**Stortz Fitting**  
Small silver plug with central pin



**Oertli blue ring to 2 pin**



**Oertli blue ring to 2 pin**



A standard cable is 3 meters long. All cables are Autoclavable, contact Altomed Limited Quality Department for Instructions.